

Set up Slack

Send your customer-service alerts to the Slack workspace your team already uses — and act on them right from Slack.

What you get

Once connected, CareMetric Breathe posts real-time, non-PHI alerts into a Slack channel your reps watch, each with a deep link into the admin console:

- **A patient reply needs a human** — an SMS/email conversation moved to the queue.
- **Voice call handoff** — the AI summarizer flagged a call (with caller sentiment).
- **SLA breach** — a conversation passed its response deadline.
- **Reminder ladder exhausted** — a patient is unresponsive after every channel.
- **NPS detractor** — a low post-delivery survey score.
- **Delivery-failure spike** — outbound messages are bouncing.
- **Operator digests** — weekly KPIs, metric alerts, stuck jobs, low stock.

And your reps can act without leaving Slack:

- **Claim** — assign the conversation to yourself.
- **Escalate** — push it to the top of the queue.
- **Snooze** — defer it for a day.
- **/pennfit queue** — see unassigned + SLA-breaching counts on demand.

Privacy

Slack messages are intentionally non-PHI: a reference id, a status, and a link — never message bodies, names, phone numbers, or clinical detail. The actual content stays in the admin console.

The easy way: one-click “Add to Slack”

If your platform administrator has registered the CareMetric Breathe Slack app, this is all you need:

- 1 Go to **Admin** → **System** → **System Configuration** and find the **Team notifications (Slack)** card.
- 2 Click **Add to Slack**.
- 3 In Slack's approval screen, choose the channel your reps watch and click **Allow**.
- 4 You're returned to System Configuration with a “Connected” banner. The bot token, workspace id, and channel are filled in for you.
- 5 Click **Send test message** to confirm a message lands in the channel.

Don't see the “Add to Slack” button working?

It says “one-click connect isn't set up by the platform yet.” Use the manual steps on the next page, or ask your platform administrator to register the Slack app (client id + secret).

The manual way: connect your own Slack app

Use this if you're bringing your own Slack app. It takes about 10 minutes. You'll create a Slack app, copy three values into System Configuration, and verify.

Step 1 — Create the Slack app

- 1 Go to api.slack.com/apps and click **Create New App** → **From scratch**.
- 2 Name it (e.g. "CareMetric Breathe"), pick your workspace, and click **Create App**.

Step 2 — Add the bot permission

- 1 In the left sidebar open **OAuth & Permissions**.
- 2 Under **Scopes** → **Bot Token Scopes**, click **Add an OAuth Scope** and add **chat:write**.
- 3 (Optional, for the slash command) also add **commands**.

Step 3 — Install & copy the bot token

- 1 At the top of **OAuth & Permissions**, click **Install to Workspace** and approve.
- 2 Copy the **Bot User OAuth Token** — it starts with `xoxb-`.
- 3 In CareMetric Breathe, paste it into **Bot token** on the Slack card and click **Save**.

Step 4 — Pick the channel & invite the bot

- 1 In Slack, open the channel your reps watch (e.g. `#cs-alerts`).
- 2 Invite the bot: type `/invite @YourAppName` in that channel.
- 3 Get the channel id: click the channel name → scroll to the bottom of the **About** tab → copy the **Channel ID** (starts with `C`).
- 4 Paste it into **Alerts channel id** in CareMetric Breathe and **Save**.

That's the minimum for alerts.

With just the bot token + channel id saved, click **Send test message**. It verifies the token, auto-detects and saves your workspace id, and posts a confirmation. The buttons + slash command need two more steps below.

Step 5 — Enable the buttons (Claim / Escalate / Snooze)

- 1 In the Slack app, open **Basic Information** → **App Credentials** and copy the **Signing Secret**.
- 2 Paste it into **Signing secret** in CareMetric Breathe and Save.
- 3 In the Slack app, open **Interactivity & Shortcuts**, turn **Interactivity** on, and set the **Request URL** to the Interactivity URL shown on the Slack card (it ends in `/resupply-api/slack/interactivity`).
- 4 Save in Slack.

Step 6 — Enable the `/pennfit` slash command (optional)

- 1 In the Slack app, open **Slash Commands** → **Create New Command**.
- 2 Command: `/pennfit`. Request URL: the Slash command URL on the Slack card (ends in `/resupply-api/slack/commands`).
- 3 Short description: "Show the CS queue." Save.
- 4 Reinstall the app if Slack prompts you to (new scope).

Step 7 — Send a test message

- 1 Back on the Slack card, click **Send test message**.
- 2 You should see "Connected to <your workspace> — test message sent" and a ■ message in your channel.
- 3 Done. Alerts will now post automatically.

Where each value comes from

Value	Where to find it
Bot token	Slack → OAuth & Permissions → Bot User OAuth Token (xoxb - ...).
Alerts channel id	Slack channel → About tab → Channel ID (C...). Invite the bot first.
Signing secret	Slack → Basic Information → App Credentials. Only needed for buttons / slash command.
Workspace (team) id	Auto-detected and saved when you click Send test message — you don't enter it.
Interactivity / Slash URLs	Shown on the Slack card; paste into the Slack app's Interactivity and Slash Command settings.

Turning features on/off

Three switches under **Admin** → **System** → **Feature Flags** control Slack (all on by default, but inert until you connect):

- **slack.notifications** — real-time CS alerts.
- **slack.digests** — the operator digests (set an #ops channel with the optional *Digests channel id* to separate them).
- **slack.interactivity** — inbound buttons + the /pennfit command.

Make “Claim” assign to the right rep

So a Claim click assigns the conversation to the person who clicked, link each rep's Slack account to their admin account:

- 1 In Slack, the rep opens their profile → **More (...)** → **Copy member ID** (starts with U).
- 2 In CareMetric Breathe, go to **Admin** → **Team**, find the rep, and paste it into the **Slack ID** field.
- 3 Now their Claim clicks assign to them. Unlinked reps get a friendly prompt to link their account.

Troubleshooting

Symptom	What it means	Fix
Test says not configured	Bot token or channel id is missing.	Add both on the Slack card, then retry.
Test says invalid_auth / token rejected	The bot token is wrong or the app was uninstalled.	Re-copy the xoxb- token; reinstall the app.
Test says channel_not_found	The channel id is wrong, or it's a private channel the bot can't see.	Re-copy the Channel ID from the About tab; use a channel the bot is in.
Test says not_in_channel	The bot isn't a member of the channel.	Run <code>/invite @YourApp</code> in that channel.
Buttons do nothing	Signing secret missing, or the Interactivity Request URL is wrong.	Paste the signing secret; set the exact Interactivity URL from the card.
<code>/pennfit</code> not found	The slash command isn't created or points at the wrong URL.	Create <code>/pennfit</code> with the Slash command URL from the card; reinstall.
Claim says "account isn't linked"	The clicking rep has no Slack ID on their Team record.	Add their Slack member id under Admin → Team.
No alerts arriving	Connected, but the <code>slack.notifications</code> flag is off.	Enable it under Admin → System → Feature Flags.

Still stuck?

Open **Admin** → **System** → **Support** and file a ticket — include what the **Send test message** button reported.

CareMetric Breathe connects to the Slack workspace your team already uses; nothing new for your reps to install.